



Low-level techniques applied in experiments looking for rare events

Grzegorz Zuzel

Max-Planck-Institut für Kernphysik, Heidelberg, Germany

Radon detection

Mass spectrometry

Germanium spectroscopy


Applications in Borexino

Applications in GERDA

Conclusions



Borexino Design




Borexino details: see the talks by G. Bellini and M. Pallavicini

1. Radon detection

Radon (^{222}Rn) and its daughters form one of the most dangerous sources of background in many experiments

- inert noble gas
- high diffusion and permeability
- belongs to the ^{238}U chain (present in any material)
- wide range of energy of emitted radiation (with the daughters)
- surface contaminations with radon daughters (heavy metals)
- broken equilibrium in the chain at ^{210}Pb level



Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino


Applications in GERDA

Conclusions



1. Radon detection

Proportional counters



- Developed for the GALLEX/GNO experiment
- Hand-made at MPIK ($\sim 1 \text{ cm}^3$ active volume)
- In case of ^{222}Rn only α -decays are detected
- 50 keV threshold
 - bcg: 0.1 – 2 cpd
 - total detection efficiency of $\sim 1.5 (0.5/\alpha)$
- **Absolute detection limit $\sim 30 \mu\text{Bq}$ (15 atoms)**

Int. Conference on Topics in Astroparticle and Underground Physics (TAUP) 2007, 11-15.09.2007 Sendai, Japan

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



1. Radon detection

Radon detection


Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



^{222}Rn in gases (N_2/Ar)

- ^{222}Rn adsorption on activated carbon
- several AC traps available (MoREx/MoRExino)
- pre-concentration from 100 – 200 m^3
- purification is possible (LTA)

^{222}Rn detection limit:
 $\sim 0.5 \mu\text{Bq}/\text{m}^3$ (STP)
[1 atom in 4 m^3]



^{222}Rn emanation

- Emanation chambers
20 l \rightarrow 50 μBq
80 l \rightarrow 80 μBq
- Glass vials
1 l \rightarrow $\sim 50 \mu\text{Bq}$

Absolute sensitivity
 $\sim 100 \mu\text{Bq}$ [50 atoms]

1. Radon detection

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions




$^{222}\text{Rn}/^{226}\text{Ra}$ in water

- ^{222}Rn extraction from 350 liters
- ^{222}Rn and ^{226}Ra measurements possible

^{222}Rn detection limit: $\sim 0.1 \text{ mBq/m}^3$

^{226}Ra detection limit: $\sim 0.8 \text{ mBq/m}^3$



^{222}Rn diffusion in thin films

- Time dependent diffusion profile registered
- D reconstructed on the base of the mathematical model

Sensitivity $\sim 10^{-13} \text{ cm}^2/\text{s}$

2. Mass spectrometry

Noble gas mass spectrometer

Radon detection


Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



VG 3600
magnetic sector field
spectrometer.

Used to investigate
noble gases in the
terrestrial and extra-
terrestrial samples.

Adopted to test the
nitrogen purity and
purification methods.

Detection limits: Ar: 10^{-9} cm³
Kr: 10^{-13} cm³



³⁹Ar: $\sim 1.4 \times 10^{-9}$ Bq/m³ N₂ (STP)
⁸⁵Kr: $\sim 1 \times 10^{-7}$ Bq/m³ N₂ (STP)

3. Germanium spectroscopy


GeMPIs at GS (3800 m w.e.)

- GeMPI I operational since 1997 (MPIK)
- GeMPI II built in 2004 (MCavern)
- GeMPI III constructed in 2007 (MPIK/LNGS)
- Worlds most sensitive spectrometers

GeMPI I:

- Crystall: 2.2 kg, $\epsilon_r = 102\%$
- Bcg. Index (0.1-2.7 MeV): 6840 cts/kg/year
- Sample chamber: 15 l

Sensitivity:
 $\sim 10 \mu\text{Bq/kg}$



Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



3. Germanium spectroscopy

Detectors at MPIK: Dario, Bruno and Corrado

Radon detection


Mass spectrometry

Germanium spectroscopy

Applications in Borexino


Applications in GERDA

Conclusions



MPIK LLL: 15 m w.e.

Sensitivity:
~1 mBq/kg



Radon detection


Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



5. Applications in Borexino - ^{226}Ra in the nylon foil

Borexino nylon foil



**1 ppt U required
(~12 $\mu\text{Bq/kg}$ for ^{226}Ra)**

$$D_{\text{dry}} = 2 \times 10^{-12} \text{ cm}^2/\text{s} \quad (d_{\text{dry}} = 7 \text{ } \mu\text{m})$$
$$D_{\text{wet}} = 1 \times 10^{-9} \text{ cm}^2/\text{s} \quad (d_{\text{wet}} = 270 \text{ } \mu\text{m})$$

$$A_{\text{dry}} = A_{\text{sf}} + 0.14 \cdot A_{\text{bulk}}$$
$$A_{\text{wet}} = A_{\text{sf}} + A_{\text{bulk}}$$

**Separation of the bulk
and surface ^{226}Ra conc.
was possible through
 ^{222}Rn emanation**

**Very sensitive technique:
($C_{\text{Ra}} \sim 10 \text{ } \mu\text{Bq/kg}$)**

**Bx IV foil: bulk $\leq 15 \text{ } \mu\text{Bq/kg}$
 surface $\leq 0.8 \text{ } \mu\text{Bq/m}^2$
 total = $(16 \pm 4) \text{ } \mu\text{Bq/kg}$ (1.3 ppt U equiv.)**

Radon detection

Mass
spectrometry

Germanium
spectroscopy

Applications in
Borexino

Applications in
GERDA

Conclusions



5. Applications in Borexino – nitrogen

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



Regular Purity Nitrogen:

- Technical 4.0 quality, not purified
- Production rate up to 100 m³/h (STP)
- ²²²Rn ~ 50 μBq/m³, Ar ~ 10ppm, Kr ~ 30 ppt

High Purity Nitrogen:

- ²²²Rn adsorption on charcoal (LTA)
- Achieved concentration < 0.3 μBq/m³
- Production rate up to 100 m³/h (STP)
- Ar and Kr not removed



LTA

LAK (Low Ar and Kr) Nitrogen:

- Spec. Ar < 0.4 ppm, Kr < 0.2 ppt
²²²Rn < 7 μBq/m³
- Purification by adsorption on different materials extensively studied (successfully!)
- Cooperation with companies on the nitrogen survey
- Tests of the nitrogen delivery chain

Nitrogen survey

Nitrogen sample	C _{Ar} [ppm]	C _{Kr} [ppt]
MESSER (4.0)	200 ± 30	1680 ± 240
Air Liquide (4.0)	11.0 ± 1.3	40 ± 5
Linde AG, (7.0)	0.031 ± 0.004	2.9 ± 0.4
SOL (6.0)	0.0063 ± 0.0006	0.04 ± 0.01
Westfalen AG (6.0)	0.00050 ± 0.00008	0.06 ± 0.02
Goal (BOREXINO)	< 0.4	< 0.2

Tests of the delivery chains



Supplier/setup	C _{Rn} [μBq/m ³]	C _{Ar} [ppm]	C _{Kr} [ppt]
Linde AG, 3-m ³ movable tank	1.2	0.018	0.06
SOL, 16-m ³ tank	8	0.012	0.02

5. Applications in Borexino – nitrogen

LAK Nitrogen tank installed at Gran Sasso

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



Int. Conference on Topics in Astroparticle and Underground Physics (TAUP) 2007, 11-15.09.2007 Sendai, Japan

6. Applications in GERDA

GERDA details: see the talk by H. Simgen

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA


Conclusions



Additional inner copper shield
detectors holders
 $U/Th \leq 16 \mu\text{Bq/kg}$
clean surface

Liquid argon
 $^{222}\text{Rn} \leq 1 \mu\text{Bq/m}^3$

Vacuum-insulated double
wall stainless steel cryostat
 $U/Th \leq 5 - 10 \text{ mBq/kg}$
clean surface



6. Applications in GERDA – argon purity

^{222}Rn in argon (GERDA goal: $C_{\text{Rn}} \leq 1 \mu\text{Bq}/\text{m}^3$)

Quality	Company	Sample size [m ³]	C_m [mBq/m ³]	C_f [mBq/m ³]
4.6	Westfalen	117	2.9 ± 0.2	>8
4.8	Air Liquide	80	0.27 ± 0.02	0.4
5.0	Westfalen	200	6.0 ± 0.1	8.4
5.0	Air Liquide (GS)	5	0.25 ± 0.04	0.3
6.0	Westfalen	104	0.11 ± 0.01	0.4

^{222}Rn in nitrogen

Quality	Company	Sample size [m ³]	C_m [$\mu\text{Bq}/\text{m}^3$]	C_f [$\mu\text{Bq}/\text{m}^3$]
4.5 – 5.0	Messer/Air Liquide	40	30 – 70	---
6.0 (7.0)	Linde	150	0.7 ± 0.2	1
6.0	SOL	100	15 ± 1	17

Argon purification required. In the gas phase requested purity achieved ($C_{\text{Rn}} \leq 0.5 \mu\text{Bq}/\text{m}^3$), adsorption for liquid phase under investigations.

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions



6. Applications in GERDA – steel screening

GERDA goal for U/Th: $\leq 5 - 10$ mBq/kg

No.	Specific activity [mBq/kg]			
	^{228}Th	^{226}Ra	^{40}K	^{60}Co
1 D	5.1 ± 1.0	2.9 ± 1.0	< 3.9	6.5 ± 0.5
2 G	< 0.27	< 0.35	< 1.1	13.0 ± 0.6
3 D	1.1 ± 0.4	< 0.84	< 3.3	15.1 ± 0.5
4 D	< 2.6	< 2.2	< 6.2	14.4 ± 1.0
5 D	< 1.1	< 1.2	< 2.8	11.6 ± 0.5
6 D	< 0.8	< 0.6	< 1.7	16.7 ± 0.4
7 G	< 0.20	< 1.3	< 2.8	45.5 ± 2.1
8 G	< 0.11	< 0.24	< 0.93	14.0 ± 0.1
9 G	< 0.41	< 0.74	< 1.1	13.8 ± 0.7
10 G	< 1.0	< 1.3	< 6.8	17.1 ± 0.7
11 G	1.5 ± 0.2	1.0 ± 0.6	< 0.81	18.3 ± 0.7

(SS 1.4571)

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA


Conclusions



6. Applications in GERDA – copper cleaning

^{222}Rn daughters on copper surface

- Screening of ^{210}Po with an alpha spectrometer
50 mm Si-detector, bcg $\sim 5 \alpha/\text{d}$ (1-10 MeV)
sensitivity $\sim 20 \text{ mBq}/\text{m}^2$ (100 mBq/kg, ^{210}Po)
- Screening of ^{210}Bi with a beta spectrometer
2×50 mm Si(Li)-detectors, bcg $\sim 0.18/0.40 \text{ cpm}$
sensitivity $\sim 10 \text{ Bq}/\text{kg}$
- Screening of ^{210}Pb (46.6 keV line) with a gamma spectrometer
25 % - n-type HPGe detector with an active and a passive shield
sensitivity $\sim 20 \text{ Bq}/\text{kg}$
- Only small samples can be handled – artificial contamination
needed: e.g. discs loaded with ^{222}Rn daughters



Copper cleaning tests

- Etching removes most of ^{210}Pb and ^{210}Bi ($> 98 \%$) but **not** ^{210}Po
- Electropolishing is more effective for all elements but proper
conditions have to be found (e.g. ^{210}Po reduction from 30 up to 200)

Etching: 1% H_2SO_4 + 3% H_2O_2 Electropolishing: 85 % H_3PO_4 + 5 % 1-butanol

Int. Conference on Topics in Astroparticle and Underground Physics (TAUP) 2007, 11-15.09.2007 Sendai, Japan

Radon detection

Mass
spectrometry

Germanium
spectroscopy

Applications in
Borexino

Applications in
GERDA

Conclusions



7. Conclusions

Introduction

Radon detection

Mass
spectrometry

Germanium
spectroscopy

Applications in
Borexino

Applications in
GERDA

Conclusions



- Low-level techniques have “natural” applications in experiments looking for rare events (low-energy neutrinos, neutrino-less double beta decay, search for dark matter/proton decay...)
- Several detectors and experimental methods were developed allowing measurements even at a single atom level
- Described experimental methods were very successfully applied in the Borexino experiment (very low background achieved) and can be adopted in other projects
 - material screening
 - purification and cleaning techniques
 - study of noble gases

6. Applications in GERDA – argon purification

^{222}Rn removal from gaseous/liquid argon

150 g- (gas phase) and 60 g-AC traps (liquid phase) used

Radon detection

Mass spectrometry

Germanium spectroscopy

Applications in Borexino

Applications in GERDA

Conclusions

Quality	Sample size [m ³]	C ₁ [mBq/m ³]	C ₂ [mBq/m ³]	Red. factor [1/kg]	Remarks
4.6	141	0.20 ± 0.02	< 0.0005	> 2700	Gas phase
4.8	80	0.27 ± 0.02	0.0007 ± 0.0003	2500	Gas phase
4.6	67	0.050 ± 0.003	0.0020 ± 0.0005	420	Liquid phase
4.6	77	0.056 ± 0.004	0.0027 ± 0.0006	370	Liquid phase
4.8	140	0.20 ± 0.01	0.005 ± 0.001	640	Liquid phase
4.8	48	0.14 ± 0.01	0.003 ± 0.001	700	Liquid phase
5.0	200	6.0 ± 0.1	0.60 ± 0.02	170	Liquid phase
6.0	104	0.11 ± 0.01	0.006 ± 0.001	305	Liquid phase

Required reduction factor for GERDA: O(500)