


PHASE I detector prototype testing...summary and status



M.Barnabé Heider, C.Cattadori, O.Chkvorets, K.Gusev, S.Schönert, M.Shirchenko


February 2007, GERDA collaboration meeting, Ringberg Castle

Outline


- Summary of 1 year measurement with Phase I detector prototype
 - Mounting
 - Cooling-warming cycles
 - Monitoring (leakage current, spectra)
 - Electronics testing
 - Refurbishment
- Long term stability measurement started

1 year testing...


- Low mass copper holder



Mounting



HV contact quality




Mechanical stability



1 year testing...

- Test stand of the LArGe facility



Nitrogen, Argon filling



Spectroscopy performance



Argon level monitoring



Infrared shielding




Radon 'free' clean bench





Radon, humidity, temperature monitoring

1 year testing...

- Electronics testing



Warm preamp




Cold preamp



HV feedthrough

Detector 'health' monitoring


- Signal to HV resistivity measurement




- Test point, picoammeter and noise level recorded



- Forward resistivity measurement while cooling





- Spectra collecting




Detector 'health' monitoring


- Latest gadget : leakage current probe (à la Canberra)
 - I-V curve of crystals



February 2007, GERDA collaboration meeting, Ringberg Castle



February 2007, GERDA collaboration meeting, Ringberg Castle



February 2007, GERDA collaboration meeting, Ringberg Castle

1 year testing...

- 43 cooling-warming cycles:
32→refurbishments, 8→refurbishments
- We know
 - Origin of leakage current: passivation layer
 - High humidity harms passivation layer
 - High cleanliness required: gloves, methanol, air, particulate, ...
- Questions
 - Handling?
 - Passivation layer?
 - Argon purity?

Conclusion

- Detector parameters measured in Canberra are reproduced in GERDA detector lab (10 pA)
- Long term stability in Argon measurement started (^{60}Co : 3.4 keV ; Pulser: 3.0 keV)
- Dedicated set up for electronics testing